



5mm Cutting Burr 30cms Shaft (Imported German Made)

5mm Cutting Burr 30cms Shaft (Imported German Made)

[Read More](#)

SKU: MCB500

Price: ₹0.00

Stock: instock

Categories: [Drills, Perforators, Trephines](#), [Endoscopic Interlaminar, Neuro & Spine, Perforators, Trephines](#)

Tags: [MCB500](#)

Product Description

5mm Cutting Burr 30cms Shaft (Imported German Made)



Autoclavable Container With Silicon Inserts & Bridge To Hold Required Set

Autoclavable Container With Silicon Inserts & Bridge To Hold Required Set

[Read More](#)

SKU: DMB404

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar, Misc Instruments, Neuro & Spine](#)

Tags: [DMB404](#), [FL109](#)

Product Description

Autoclavable Container With Silicon Inserts & Bridge To Hold Required Set



Curette Fibre Handle 3mm Cup Working Length 30cms Curved

Curette Fibre Handle 3mm Cup Working Length 30cms Curved

[Read More](#)

SKU: DMB523

Price: ₹0.00

Stock: instock

Categories: [Curettes](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB523](#), [FL553](#)

Product Description

Curette Fibre Handle 3mm Cup Working Length 30cms Curved



Curette Fibre Handle 3mm Cup Working Length 30cms Straight

Curette Fibre Handle 3mm Cup Working Length 30cms Straight

[Read More](#)

SKU: DMB521

Price: ₹0.00

Stock: instock

Categories: [Curettes](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB521](#)

Product Description

Curette Fibre Handle 3mm Cup Working Length 30cms Straight



Dilator Single Bore 10mm 23cms With Crimped Handle Grip

Dilator Single Bore 10mm 23cms With Crimped Handle Grip

[Read More](#)

SKU: DMB501

Price: ₹0.00

Stock: instock

Categories: [Dilator](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB501](#), [FL505](#)

Product Description

Dilator Single Bore 10mm 23cms With Crimped Handle Grip



Disc Rongeur Tubular Wl 30cms
Shaft 4mm 3mm X 10mm Jaw
Serr Straight

Disc Rongeur Tubular Wl 30cms Shaft 4mm 3mm X 10mm Jaw Serr
Straight

[Read More](#)

SKU: DMB541

Price: ₹0.00

Stock: instock

Categories: [Disc Ronguers](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB541](#), [FL583](#)

Product Description

Disc Rongeur Tubular Wl 30cms Shaft 4mm 3mm X 10mm Jaw Serr Straight



Disc Rongeur Tubular Wl 30cms
Shaft 4mm 3mm X 10mm Jaw
Serr Up Curved

Disc Rongeur Tubular Wl 30cms Shaft 4mm 3mm X 10mm Jaw Serr Up
Curved

[Read More](#)

SKU: DMB543

Price: ₹0.00

Stock: instock

Categories: [Disc Ronguers](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB543](#), [FL585](#)

Product Description

Disc Rongeur Tubular Wl 30cms Shaft 4mm 3mm X 10mm Jaw Serr Up Curved



Disc Rongeur Tubular Wl 30cms Shaft 4mm 4mm X 10mm Jaw Serr Up Curved

Disc Rongeur Tubular Wl 30cms Shaft 4mm 4mm X 10mm Jaw Serr Up Curved

[Read More](#)

SKU: DMB545

Price: ₹0.00

Stock: instock

Categories: [Disc Ronguers](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB545](#)

Product Description

Disc Rongeur Tubular Wl 30cms Shaft 4mm 4mm X 10mm Jaw Serr Up Curved



DissectorCvd Fibre Handle 3mm Tip Working Length 30cms

DissectorCvd Fibre Handle 3mm Tip Working Length 30cms

[Read More](#)

SKU: DMB517

Price: ₹0.00

Stock: instock

Categories: [Dissectors](#), [Raspatory](#), [Elevators etc](#), [Elevators etc](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#), [Raspatory](#)

Tags: [DMB517](#), [FL526](#)

Product Description

DissectorCvd Fibre Handle 3mm Tip Working Length 30cms



Endoscope Interlaminar 10mm 12.5cms15 Degree / 25 Degree Imported

Endoscope Interlaminar 10mm 12.5cms15 Degree / 25 Degree Imported

[Read More](#)

SKU: ESC200

Price: ₹0.00

Stock: instock

Categories: [Endoscope](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)
Tags: [ESC200](#)

Product Description

Endoscope Interlaminar 10mm 12.5cms15 Degree / 25 Degree Imported



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm

[Read More](#)

SKU: DMB529
Price: ₹0.00
Stock: instock
Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)
Tags: [DMB529](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 2mm



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm

[Read More](#)

SKU: DMB531
Price: ₹0.00
Stock: instock
Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)
Tags: [DMB531](#), [FL573](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 3mm



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm

[Read More](#)

SKU: DMB533

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB533](#), [FL575](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 4mm



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm

[Read More](#)

SKU: DMB535

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB535](#), [FL576](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 40 Deg Up 5mm



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up 4mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up 4mm

[Read More](#)

SKU: DMB539

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB539](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up 4mm



Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm

[Read More](#)

SKU: DMB537

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Kerrison Ronguers](#), [Neuro & Spine](#)

Tags: [DMB537](#), [FL577](#)

Product Description

Kerrison Rongeur Tubular Wl 30cms Shaft 5mm 90 Deg Up3mm



Mallet Cottle 250 G

Mallet Cottle 250 G

[Read More](#)

SKU: DMB505

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Mallet](#), [Neuro & Spine](#)

Tags: [DMB505](#), [FL507](#)

Product Description

Mallet Cottle 250 G



Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd

Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd

[Read More](#)

SKU: DMB527

Price: ₹0.00

Stock: instock

Categories: [Chisels](#), [Gouges](#), [Neuro & Spine](#), [Endoscopic Interlaminar](#), [Osteotome](#)

Tags: [DMB527](#), [FL556](#)

Product Description

Osteotome 4mm Tip Working Length 30cms Fibre Handle Cvd



Osteotome 4mm Tip Working Length 30cms Fibre Handle St

Osteotome 4mm Tip Working Length 30cms Fibre Handle St

[Read More](#)

SKU: DMB525

Price: ₹0.00

Stock: instock

Categories: [Chisels](#), [Endoscopic Interlaminar](#), [Gouges](#), [Neuro & Spine](#), [Osteotome](#)

Tags: [DMB525](#), [FL555](#)

Product Description

Osteotome 4mm Tip Working Length 30cms Fibre Handle St



Right Angle Ball Probe Fibre
Handle 3mm Working Length
30cms

Right Angle Ball Probe Fibre Handle 3mm Working Length 30cms

[Read More](#)

SKU: DMB519
Price: ₹0.00
Stock: instock
Categories: [Endoscopic Interlaminar](#), [Hooks & Probes](#), [Neuro & Spine](#)
Tags: [DMB519](#), [FL527](#)

Product Description

Right Angle Ball Probe Fibre Handle 3mm Working Length 30cms



T Handle Cannula 11.2mm Od
10.20mm Id Working Length
12cms Elevator Tip 10mm

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms
Elevator Tip 10mm

[Read More](#)

SKU: DMB510
Price: ₹0.00
Stock: instock
Categories: [Endoscopic Interlaminar](#), [Misc Instruments](#), [Neuro & Spine](#)
Tags: [DMB510](#), [FL525](#)

Product Description

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Elevator Tip 10mm



T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Elevator Tip 5mm

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms
Elevator Tip 5mm

[Read More](#)

SKU: DMB509

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB509](#), [FL523](#)

Product Description

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Elevator Tip 5mm



T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Zero Degree

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Zero
Degree

[Read More](#)

SKU: DMB507

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB507](#)

Product Description

T Handle Cannula 11.2mm Od 10.20mm Id Working Length 12cms Zero Degree



Tapper AluminiumId 11mm

Tapper AluminiumId 11mm

[Read More](#)

SKU: DMB503

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB503](#), [FL506](#)

Product Description

Tapper AluminiumId 11mm



Thrucut Forceps Tubular Wl 30cms Shaft 3.5mm 3mm X 8mm 45 Degree Angled Up

Thrucut Forceps Tubular Wl 30cms Shaft 3.5mm 3mm X 8mm 45 Degree Angled Up

[Read More](#)

SKU: DMB549

Price: ₹0.00

Stock: instock

Categories: [Bone Cutters](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB549](#)

Product Description

Thrucut Forceps Tubular Wl 30cms Shaft 3.5mm 3mm X 8mm 45 Degree Angled Up



Thrucut Forceps Tubular Wl 30cms Shaft 4mm 3mm X 10mm Upcurved

Thrucut Forceps Tubular Wl 30cms Shaft 4mm 3mm X 10mm Upcurved

[Read More](#)

SKU: DMB547

Price: ₹0.00

Stock: instock

Categories: [Bone Cutters](#), [Endoscopic Interlaminar](#), [Neuro & Spine](#)

Tags: [DMB547](#)

Product Description

Thru-cut Forceps Tubular W/ 30cms Shaft 4mm 3mm X 10mm Upcurved



Tubular Nerve Root Retractor And Gouge With Handle W/ 30cms 4mm Od 3.5mm Id

Tubular Nerve Root Retractor And Gouge With Handle W/ 30cms 4mm
Od 3.5mm Id

[Read More](#)

SKU: DMB528

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Interlaminar](#), [Hand Held Retractor](#),
[Spatulas](#), [Neuro & Spine](#), [Spatulas](#)

Tags: [DMB528](#), [FL557](#)

Product Description

Tubular Nerve Root Retractor And Gouge With Handle W/ 30cms 4mm Od 3.5mm Id