



Autoclavable Container With Silicon Inserts To Hold Transilliac Set

Autoclavable Container With Silicon Inserts To Hold Transilliac Set

[Read More](#)

SKU: DMB402

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Transilliac](#), [Misc Instruments](#), [Neuro & Spine](#)

Tags: [DMB402](#), [FL106](#)

Product Description

Autoclavable Container With Silicon Inserts To Hold Transilliac Set



Probe With Cannula & Perforated Round Dismantable Handle

Probe With Cannula & Perforated Round Dismantable Handle

[Read More](#)

SKU: DMB292

Price: ₹0.00

Stock: instock

Categories: [Endoscopic Transilliac](#), [Hooks & Probes](#), [Neuro & Spine](#)

Tags: [DMB292](#)

Product Description

Probe With Cannula & Perforated Round Dismantable Handle



Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm

[Read More](#)

SKU: DMB297

Price: ₹0.00

Stock: instock

Categories: [Perforators, Trephines, Neuro & Spine, Endoscopic Transilliac, Drills, Perforators, Trephines](#)

Tags: [DMB297](#)

Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm



Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm

[Read More](#)

SKU: DMB293

Price: ₹0.00

Stock: instock

Categories: [Drills, Perforators, Trephines, Endoscopic Transilliac, Neuro & Spine, Perforators, Trephines](#)

Tags: [DMB293](#)

Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 6mm



Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm

[Read More](#)

SKU: DMB295

Price: ₹0.00

Stock: instock

Categories: [Drills, Perforators, Trephines, Endoscopic Transilliac, Neuro & Spine, Perforators, Trephines](#)

Tags: [DMB295](#)

Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm
