




DeBakey Bull Dog Clamps 10 Cm-4"

DeBakey Bull Dog Clamps 10 Cm-4"

[Read More](#)

SKU: VBD070

Price: ₹0.00


Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.02.10](#)

Product Description

DeBakey Bull Dog Clamps 10 Cm-4"



DeBakey Bull Dog Clamps 10.5 Cm-4 1/4"

DeBakey Bull Dog Clamps 10.5 Cm-4 1/4"

[Read More](#)

SKU: VBD050

Price: ₹0.00


Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.01.10](#)

Product Description

DeBakey Bull Dog Clamps 10.5 Cm-4 1/4"



DeBakey Bull Dog Clamps 11.5 Cm- 4 3/4"

DeBakey Bull Dog Clamps 11.5 Cm- 4 3/4"

[Read More](#)

SKU: VBD075

Price: ₹0.00

Stock: instock

Categories: [Vascular](#), [Bull Dog Clamps](#)

Tags: [33.02.12](#)

Product Description

Debakey Bull Dog Clamps 11.5 Cm- 4 3/4"



Debakey Bull Dog Clamps 12Cm - 4 3/4"

Debakey Bull Dog Clamps 12Cm - 4 3/4"

[Read More](#)

SKU: VBD055

Price: ₹0.00


Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.01.12](#)

Product Description

Debakey Bull Dog Clamps 12Cm - 4 3/4"



Debakey Bull Dog Clamps 7 Cm- 2 3/4"

Debakey Bull Dog Clamps 7 Cm- 2 3/4"

[Read More](#)

SKU: VBD060

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.02.07](#)

Product Description

Debakey Bull Dog Clamps 7 Cm- 2 3/4"



DeBakey Bull Dog Clamps 7.5 Cm-3"

DeBakey Bull Dog Clamps 7.5 Cm-3"

[Read More](#)

SKU: VBD040

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.01.07](#)

Product Description

DeBakey Bull Dog Clamps 7.5 Cm-3"



DeBakey Bull Dog Clamps 8 Cm- 3 1/4"

DeBakey Bull Dog Clamps 8 Cm- 3 1/4"

[Read More](#)

SKU: VBD065

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.02.08](#)

Product Description

DeBakey Bull Dog Clamps 8 Cm- 3 1/4"



DeBakey Bull Dog Clamps 8.5 Cm- 3 1/2"

DeBakey Bull Dog Clamps 8.5 Cm- 3 1/2"

[Read More](#)

SKU: VBD045

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [33.01.08](#)

Product Description

Debakey Bull Dog Clamps 8.5 Cm- 3 1/2"



Debakey Mini Bull Dog Clamps
Ang 55mm-2 1/4"

Debakey Mini Bull Dog Clamps Ang 55mm-2 1/4"

[Read More](#)

SKU: VBD025

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD025](#)

Product Description

Debakey Mini Bull Dog Clamps Ang 55mm-2 1/4"



Debakey Mini Bull Dog Clamps
Cvd 55mm-2 1/4"

Debakey Mini Bull Dog Clamps Cvd 55mm-2 1/4"

[Read More](#)

SKU: VBD030

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD030](#)

Product Description

Debakey Mini Bull Dog Clamps Cvd 55mm-2 1/4"



Debakey Mini Bull Dog Clamps Str 55mm-2 1/4"

Debakey Mini Bull Dog Clamps Str 55mm-2 1/4"

[Read More](#)

SKU: VBD020

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD020](#)

Product Description

Debakey Mini Bull Dog Clamps Str 55mm-2 1/4"



Diethrich Bull Dog / Spikey Clamps Ang 2 1/4"

Diethrich Bull Dog / Spikey Clamps Ang 2 1/4"

[Read More](#)

SKU: VBD035

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD035](#)

Product Description

Diethrich Bull Dog / Spikey Clamps Ang 2 1/4"



Glover Bull Dog Clamps Ang 35mm

Glover Bull Dog Clamps Ang 35mm

[Read More](#)

SKU: VBD010

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD010](#)

Product Description

Glover Bull Dog Clamps Ang 35mm



Glover Bull Dog Clamps Cvd 35mm

Glover Bull Dog Clamps Cvd 35mm

[Read More](#)

SKU: VBD015

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD015](#)

Product Description

Glover Bull Dog Clamps Cvd 35mm



Glover Bull Dog Clamps Str 35mm

Glover Bull Dog Clamps Str 35mm

[Read More](#)

SKU: VBD005

Price: ₹0.00

Stock: instock

Categories: [Bull Dog Clamps](#), [Vascular](#)

Tags: [VBD005](#)

Product Description

Glover Bull Dog Clamps Str 35mm