



Abdominal Trochar Blunt Tip Channeled Sheeth

Abdominal Trochar Blunt Tip Channeled Sheeth

[Read More](#)

SKU: NTR005

Price: ₹0.00

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Trocar](#)

Tags: [NTR005](#)

Product Description

Abdominal Trochar Blunt Tip Channeled Sheeth



22.01.22 2.2 mm
22.01.30 3.0 mm
22.01.40 4.0 mm
22.01.50 5.0 mm
22.01.55 5.5 mm

Nelson

Nelson Trocar 20 Cm-8" 2.2 Mm

Nelson Trocar 20 Cm-8" 2.2 Mm

[Read More](#)

SKU: UTR005

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.01.22](#)

Product Description

Nelson Trocar 20 Cm-8" 2.2 Mm



22.01.22 2.2 mm
22.01.30 3.0 mm
22.01.40 4.0 mm
22.01.50 5.0 mm
22.01.55 5.5 mm

Nelson

Nelson Trocar 20 Cm-8" 3. Mm

Nelson Trocar 20 Cm-8" 3. Mm

[Read More](#)

SKU: UTR010

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.01.30](#)

Product Description

Nelson Trocar 20 Cm-8" 3. Mm



Nelson Trocar 20 Cm-8" 4.0 Mm

Nelson Trocar 20 Cm-8" 4.0 Mm

[Read More](#)

SKU: UTR015
Price: ₹0.00
Stock: instock
Categories: [Trocar](#), [Urology](#)
Tags: [22.01.40](#)

Product Description

Nelson Trocar 20 Cm-8" 4.0 Mm



Nelson Trocar 20 Cm-8" 5.0 Mm

Nelson Trocar 20 Cm-8" 5.0 Mm

[Read More](#)

SKU: UTR020
Price: ₹0.00
Stock: instock
Categories: [Urology](#), [Trocar](#)
Tags: [22.01.50](#)

Product Description

Nelson Trocar 20 Cm-8" 5.0 Mm



22.01.22 2.2 mm
22.01.30 3.0 mm
22.01.40 4.0 mm
22.01.50 5.0 mm
22.01.55 5.5 mm

Nelson

Nelson Trocar 20 Cm-8" 5.5 Mm

Nelson Trocar 20 Cm-8" 5.5 Mm

[Read More](#)

SKU: UTR025

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.01.55](#)

Product Description

Nelson Trocar 20 Cm-8" 5.5 Mm



22.02.30 3.0 mm
22.02.36 3.6 mm
22.02.43 4.3 mm
22.02.50 5.0 mm

Ochsner

Ochsner Trocar 20 Cm-8" 3.0 Mm

Ochsner Trocar 20 Cm-8" 3.0 Mm

[Read More](#)

SKU: UTR030

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.02.30](#)

Product Description

Ochsner Trocar 20 Cm-8" 3.0 Mm



22.02.30 3.0 mm
22.02.36 3.6 mm
22.02.43 4.3 mm
22.02.50 5.0 mm

Ochsner

Ochsner Trocar 20 Cm-8" 3.6 Mm

Ochsner Trocar 20 Cm-8" 3.6 Mm

[Read More](#)

SKU: UTR035

Price: ₹0.00


Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.02.36](#)

Product Description

Ochsner Trocar 20 Cm-8" 3.6 Mm



Ochsner Trocar 20 Cm-8" 4.3 Mm


Ochsner Trocar 20 Cm-8" 4.3 Mm

[Read More](#)

SKU: UTR040
Price: ₹0.00
Stock: instock
Categories: [Trocar](#), [Urology](#)
Tags: [22.02.43](#)

Product Description

Ochsner Trocar 20 Cm-8" 4.3 Mm



Ochsner Trocar 20 Cm-8" 5.0 Mm

Ochsner Trocar 20 Cm-8" 5.0 Mm

[Read More](#)

SKU: UTR045
Price: ₹0.00
Stock: instock
Categories: [Trocar](#), [Urology](#)
Tags: [22.02.50](#)

Product Description

Ochsner Trocar 20 Cm-8" 5.0 Mm

Trocar Button Model

Trocar Button Model



[Read More](#)

SKU: UTR050

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [UTR050](#)

Product Description

Trocar Button Model
