



# Abdominal Trochar Blunt Tip Channeled Sheeth

Abdominal Trochar Blunt Tip Channeled Sheeth

[Read More](#)

**SKU:** NTR005

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Neuro & Spine](#), [Open Surgery](#), [Trocar](#)

**Tags:** [NTR005](#)

## Product Description

Abdominal Trochar Blunt Tip Channeled Sheeth



22.01.22 2.2 mm  
22.01.30 3.0 mm  
22.01.40 4.0 mm  
22.01.50 5.0 mm  
22.01.55 5.5 mm

Nelson

# Nelson Trocar 20 Cm-8" 2.2 Mm

Nelson Trocar 20 Cm-8" 2.2 Mm

[Read More](#)

**SKU:** UTR005

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [22.01.22](#)

## Product Description

Nelson Trocar 20 Cm-8" 2.2 Mm



22.01.22 2.2 mm  
22.01.30 3.0 mm  
22.01.40 4.0 mm  
22.01.50 5.0 mm  
22.01.55 5.5 mm

Nelson

# Nelson Trocar 20 Cm-8" 3. Mm

Nelson Trocar 20 Cm-8" 3. Mm

[Read More](#)

**SKU:** UTR010

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [22.01.30](#)

Product Description

Nelson Trocar 20 Cm-8" 3. Mm



Nelson Trocar 20 Cm-8" 4.0 Mm

Nelson Trocar 20 Cm-8" 4.0 Mm

[Read More](#)

SKU: UTR015

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.01.40](#)

Product Description

Nelson Trocar 20 Cm-8" 4.0 Mm



Nelson Trocar 20 Cm-8" 5.0 Mm

Nelson Trocar 20 Cm-8" 5.0 Mm

[Read More](#)

SKU: UTR020

Price: ₹0.00

Stock: instock

Categories: [Trocar](#), [Urology](#)

Tags: [22.01.50](#)

Product Description

Nelson Trocar 20 Cm-8" 5.0 Mm



22.01.22 2.2 mm  
22.01.30 3.0 mm  
22.01.40 4.0 mm  
22.01.50 5.0 mm  
22.01.55 5.5 mm

Nelson

## Nelson Trocar 20 Cm-8" 5.5 Mm

Nelson Trocar 20 Cm-8" 5.5 Mm

[Read More](#)

**SKU:** UTR025

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [22.01.55](#)

### Product Description

Nelson Trocar 20 Cm-8" 5.5 Mm



22.02.30 3.0 mm  
22.02.36 3.6 mm  
22.02.43 4.3 mm  
22.02.50 5.0 mm

Ochsner

## Ochsner Trocar 20 Cm-8" 3.0 Mm

Ochsner Trocar 20 Cm-8" 3.0 Mm

[Read More](#)

**SKU:** UTR030

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [22.02.30](#)

### Product Description

Ochsner Trocar 20 Cm-8" 3.0 Mm



22.02.30 3.0 mm  
22.02.36 3.6 mm  
22.02.43 4.3 mm  
22.02.50 5.0 mm

Ochsner

## Ochsner Trocar 20 Cm-8" 3.6 Mm

Ochsner Trocar 20 Cm-8" 3.6 Mm

[Read More](#)

**SKU:** UTR035

**Price:** ₹0.00


**Stock:** instock

**Categories:** [Urology](#), [Trocar](#)

**Tags:** [22.02.36](#)

Product Description

Ochsner Trocar 20 Cm-8" 3.6 Mm



Ochsner Trocar 20 Cm-8" 4.3 Mm


Ochsner Trocar 20 Cm-8" 4.3 Mm

[Read More](#)

**SKU:** UTR040  
**Price:** ₹0.00  
**Stock:** instock  
**Categories:** [Trocar](#), [Urology](#)  
**Tags:** [22.02.43](#)

Product Description

Ochsner Trocar 20 Cm-8" 4.3 Mm



Ochsner Trocar 20 Cm-8" 5.0 Mm

Ochsner Trocar 20 Cm-8" 5.0 Mm

[Read More](#)

**SKU:** UTR045  
**Price:** ₹0.00  
**Stock:** instock  
**Categories:** [Trocar](#), [Urology](#)  
**Tags:** [22.02.50](#)

Product Description

Ochsner Trocar 20 Cm-8" 5.0 Mm

# Trocar Button Model

Trocar Button Model



[Read More](#)

**SKU:** UTR050

**Price:** ₹0.00

**Stock:** instock

**Categories:** [Trocar](#), [Urology](#)

**Tags:** [UTR050](#)

## Product Description

Trocar Button Model