



Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 3mm

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working
Length 15cm 3mm

[Read More](#)

SKU: NSU150

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU150](#)

Product Description

Suction Handle Tapered 3 hole Variable Non Trauma ANG Working Length 15cm 3mm