



Dilator 6.8mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.8mm EU Scope

Dilator 6.8mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.8mm EU Scope

[Read More](#)

SKU: DMB305A

Price:

Stock: instock

Categories: [Dilator](#), [Endoscopic Transforaminal](#), [Neuro & Spine](#)

Tags: [DMB305A](#)

Product Description

Dilator 6.8mm 24 Cms Id 1.5mm Single Bore Taper Tip 6.8mm EU Scope