



# Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm

[Read More](#)

**SKU:** DMB297

**Price:**

**Stock:** instock

**Categories:** [Drills, Perforators, Trephines](#), [Endoscopic Transilliac, Neuro & Spine, Perforators, Trephines](#)

**Tags:** [DMB297](#)

## Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 10mm