



Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm

[Read More](#)

SKU: DMB295

Price:

Stock: instock

Categories: [Drills, Perforators, Trephines](#), [Endoscopic Transilliac, Neuro & Spine, Perforators, Trephines](#)

Tags: [DMB295](#)

Product Description

Reamer 2mm Cannulated Side Cutting Blunt Tip 23 Cms 8mm