



Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Length 15cm 1mm

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working
Length 15cm 1mm

[Read More](#)

SKU: NSU205

Price:

Stock: instock

Categories: [Neuro & Spine](#), [Open Surgery](#), [Suction](#)

Tags: [NSU205](#)

Product Description

Suction Handle Tapered Key Hole Variable Non Trauma ANG Working Length 15cm 1mm